

ALS beamline 12.2.2




Martin Kunz


Outline:

- Beamline layout.
- Experimental station.
- Diffraction experiments: Resistive heating.
- X-ray absorption imaging (Walker/Walter).
- Viscosity measurements (Grocholski/Jeanloz).
- Radial diffraction (Merkel/Lowell/Wenk).
- Development: EXAFS.
- Development: In situ laser heated radial diffraction.
- Development: Upgrade end-station 1.


Beamline 12.2.2: Layout



2006 Calipso Review: Beamline Layout




Beamline 12.2.2: Monochromator




2006 Calipso Review: Beamline Layout

Beamline 12.2.2: Flux



2006 Calipso Review: Experimental Stations.

End station 1:



Beam position monitor 1

Filters

Shutter and 2nd bpm

4 pairs of Slits (also virtual source for K-B-mirrors)

2006 Calipso Review: Experimental Stations.



K-B-mirrors:

Virtual source is re-imaged on sample position of ES-2.

Beam size $\sim 10 \mu\text{m} \times 10 \mu\text{m}$ (theoretical).

No tails thanks to off-set slits around virtual source.

2 Si mirrors with Rh ($\sim 40 \text{ \AA}$) on Pt ($\sim 250 \text{ \AA}$) on Cr ($\sim 50 \text{ \AA}$)






CdWO₃ scintillator
~132.86mm from KB box.

KB pitches & translations
optimized earlier with
different imaging system.

Field of view is 0.920 mm
and image has 512x512
pixels => 1.80
microns/pixel


17keV



100 μm pinhole mounted immediately in front of carbon mirrors removes scatter-tails from KB-mirrors.

=> Virtually no gasket signal even for very small gasket-holes.

2006 Calipso Review: Experimental Stations.




2006 Calipso Review: Experimental Stations.





Virtual Source: 20 (h) x 60 (v) mm². Focused spot: 13 x 13 mm².


KB-Flux comparison




Resistive heating:




Resistive heating:



Resistive heating:




Zircon (ZrSiO_4):





X-ray imaging as a complementary tool.


- Measuring viscosity of liquid with a rolling sphere (Jeanloz group).
- Measuring shift of eutectic / peritectic with P (Walker and Walter).



X-ray-absorption imaging experiment using the Apogee camera



Observed absorption in FeS







50 kbars





Melting shows erosion








Observed absorption in FeO




2006 Calipso Review: Rolling Sphere Imaging



2006 Calipso Review: Imaging.






Radial x-ray diffraction geometry (Sebastien Merkel, UCB)




- Cylindrical symmetry
- Diffraction in radial geometry
- Diffraction patterns as a function of azimuthal angle δ
- One image for each deformation / pressure step
- Detector: image plate or CCD
- Data collection time: 30s (CCD) to 10mn (image plate)

2006 Calipso Review: Radial Diffraction



Merkel and Yagi, RSI (2005)


2006 Calipso Review: Radial Diffraction.



Diffraction




Preferred orientations

Texture analysis





Polycrystalline plasticity models

Orientation distribution function





Deformation mechanisms

Elasticity and stress



EXAFS





Monochromator rocking curve feedback system to be completed by end of January.

Remote controlled DAC for in-situ
laser heating of radial diffraction.



Installation of single crystal goniometer on ES-1:



Specifications:

- 4-circle goniometer with full Eulerian cradle.
- small and sturdy design to fit in available space and allowing mounting of DAC's.
- total sphere of confusion ≤ 0.01 mm radius with DAC mounted.
- vertical diffraction plane to account for horizontal polarization of synchrotron beam.
- Bruker CCD detector controlled through WinView. Possibility to add / exchange with single point detector.
- omega range $> \pm 40^\circ$.
- 2θ range $> \pm 35^\circ$.
- Controller & driver: Newport XPS.
- flexibility wrt application-software.

Further points to consider:

- Potential of technique can only be fully exploited, if gas loading (He or Ne) is available!
- Software needs to be flexible AND user-friendly. Ideally experiment can be run by user-preferred software.
- capability to expand to other ancillary equipment (heater, heatable DAC) .
- beam-size.

Resources needed.

- **Finances** for purchase of goniometer: available ~ 70 k\$ from NSF- MRI (Nicol - Dera).
- **Mechanical** support for installation (Titus / Celestre / Morrison - ALS).
- **Electronic** support for installation (Bell - ALS).
- **Controller** software (EPICS) support (Mark Rivers, GSECars).
- **Detector** interfacing: Ed Domning / Rich Celestre - ALS, Mark Rivers - GSECars.
- **Application** software (IDL based) support (Przemek Dera, GSECars. R. Angel VPI (?)).
- **Beamline** support (Kunz - COMPRES / ALS).

2006 Calipso Review: Developments



Spring shut down

Nov. 06

May 07

