

## Other high-pressure facilities at the ALS:






**Martin Kunz**

## Outline:

- 11.3.1
- 4.3.1
- High-pressure lab.

- ? range 0.6 to 2.0 Å
- Focus 200 x 90  $\mu\text{m}$ , collimated down to 100 x 100  $\mu\text{m}$
- Detector: Bruker AXS Pt200 CCD
- Flux:  $10^{10}$  photons/sec at sample location.





*Si<111> Channel-Cut Monochromator*



Internal View



External View

PHYSICAL REVIEW B 74, 115405 (2006)

## Compressibility of zinc sulfide nanoparticles

B. Gilbert,<sup>1</sup> H. Zhang,<sup>2</sup> B. Chen,<sup>2</sup> M. Kunz,<sup>3</sup> F. Huang,<sup>4</sup> and J. F. Banfield<sup>2</sup>

<sup>1</sup>*Earth Sciences Division, Lawrence Berkeley National Laboratory, Berkeley, California 94720, USA*



<sup>2</sup>*Department of Earth & Planetary Sciences, University of California–Berkeley, Berkeley, California 94720, USA*


<sup>3</sup>*Advanced Light Source, Lawrence Berkeley National Laboratory, Berkeley, California 94720, USA*


<sup>4</sup>*Fujian Institute of Research on the Structure of Matter, Chinese Academy of Sciences, Fuzhou,  
Fujian 350002, People's Republic of China*

(Received 17 March 2006; revised manuscript received 18 July 2006; published 12 September 2006)














## 2006 Calipso Review: 4.3.1




## 2006 Calipso Review: 4.3.1



Gas flushing port



Talc




10Å phase

Pawley, Clark et al.



## 2006 Calipso Review: 4.3.1



**10-Å phase:**  
 **$0.9 \text{ cm}^{-1}\text{GPa}^{-1}$**   
**(2 - 10 GPa)**

**0 - 2 GPa: no change**

### Summary Talc:


- Positive P-dependence = shortening and strengthening of O-H bond, due to compression of 2:1 layer.
- Positive P-dependence means no hydrogen-bonding, as compression would cause strengthening of hydrogen bonds and weakening of O-H bonds.
- Small frequency shift because most compression taken up by weakly bonded interlayer.

### Summary 10 Å phase:

- No shift up to 2 GPa  $\Rightarrow$  no compression of 2:1 layer. All compression must be of interlayer.
- Above 2 GPa, P-dependence is similar to talc  $\Rightarrow$  similar compressional behaviour of 2:1 layer.
- No hydrogen bonding to interlayer H<sub>2</sub>O.




## Sample preparation:



## Mechanical driller:



### Electro erosion driller:



### Cryogenic loading (Ar, N<sub>2</sub>):



## Brillouin spectroscopy:

