

Micro X-Ray Absorption Spectroscopy for Environmental Science at Beamline 10.3.2


Generic Procedure

1. X-Ray Fluorescence Mapping (XRF)

Raster the X-ray beam over the sample.


Map intensity of X-ray fluorescence from different elements to locate interesting areas.

Scan sample under beam, record X-ray fluorescence as function of position.




2. X-Ray Absorption Spectroscopy (XAS)

Pick a spot and look at it with EXAFS, XANES



3. New capability: XRD

Identify mineral with which metal is associated by detecting X-rays diffracted from matrix.



Results:
Metal location and association.
Identities of dominant minerals.
Identity of sorbent minerals.
Chemical state of metal.