


# **DIFFRACTION GRATING MONOCHROMATORS (PGMs)**

**Malcolm R. Howells**

**Advanced Light Source**



# PLANE GRATING SYSTEMS: GENERAL



Eliminating  $\alpha$  between

$$\frac{m\lambda}{d_0} = \sin \alpha + \sin \beta \quad \text{and} \quad c_{ff} = \frac{\cos \beta}{\cos \alpha}$$

$$\sin \beta = \frac{\frac{m\lambda}{d} - \sqrt{\frac{m\lambda}{d}^2 \frac{1}{c_{ff}^2} + 1 - \frac{1}{c_{ff}^2}}}{1 - \frac{1}{c_{ff}^2}}$$

Focus condition for a spherical grating

$$\frac{\cos^2 \alpha}{r} - \frac{\cos \alpha}{R} + \frac{\cos^2 \beta}{r} - \frac{\cos \beta}{R} = 0 \quad \text{let } R \quad r = -r \frac{\cos^2 \beta}{\cos^2 \alpha}$$

Note that if  $r =$  then  $r =$  irrespective of  $c_{ff}$  i.e. the grating then does no focusing

Using  $c_{ff} = \frac{\cos \beta}{\cos \alpha}$  and recalling that  $M = -\frac{r}{R} \frac{\cos \alpha}{\cos \beta}$


$$M = c_{ff} \quad (\text{of the grating alone})$$



# COLLIMATED LIGHT SX700 (CLSX700)



Kunz 1968  
Petersen 1982  
Follath 1997




- This is the current and deservedly the most popular of many versions (used in ALS MES project)
- The original (Petersen) held fixed focus by maintaining  $c_{ff} = \cos \beta / \cos \alpha = 2.25$  and introduced the still-essential mechanism for positioning the plane mirror (next slide)
- CLSX700 has constant focus position *irrespective* of the value of  $c_{ff}$  which therefore becomes a user-controlled variable
- This can be used to track the efficiency maximum, suppress high orders or maximize resolution
- Sagittal focusing mirrors reduce sensitivity to manufacturing errors (Jark 1992) see next slide, other types were paraboloid (Kunz 1968), ellipsoid (Petersen 1982), elliptical cylinder (Nyholm 1985) or sphere (Padmore 1989)



# FORGIVENESS FACTOR

## TANGENTIAL-PLANE DIAGRAM



Requirement:  $\Delta s_T \quad \frac{1}{2} s_T \quad 2r \Delta_T \quad \frac{1}{2} \frac{r}{r} s_T$

$$\Delta_T \quad \frac{1}{4} \frac{s_T}{r}$$

$\Delta s_S \quad \frac{1}{2} s_S \quad 2r \theta \Delta_S \quad \frac{1}{2} \frac{r}{r} s_S$

$$\Delta_S \quad \frac{1}{4} \frac{s_S}{r\theta}$$

Reasons for the factor 2's:


1. Ray deviation is twice that of the mirror
2. Quadratic sum of  $\sigma$  and  $\sigma/2$  is  $1.1\sigma$  implying a 10% error which we declare acceptable



# THE ZEISS MIRROR MECHANISM



- Reimer 1991
- Pimpale 1991



- This mechanism, patented by BESSY and Zeiss, is a major reason for the success of this design
- The idea is to find a point about which to rotate the mirror such that the central ray always hits the grating pole - this can be done with only a few microns error
- The easy and less accurate solution is continue the mirror surfaces at maximum and minimum angle to intersect the y axis and take the average of the two points - this will be near A ( $0, D/2$ ) for small angles
- More accurate solutions (Pimpale 1991) are close to B ( $0, 3D/2$ ) - best solutions are good within  $10\mu\text{m}$  if  $D = 8\text{ cm}$  and the grazing angle is less than  $10^\circ$  - B itself is a good solution to within  $10\mu\text{m}$  if  $D = 3\text{ cm}$  and the grazing angle is less than  $10^\circ$



# SOME WORKING CURVES FOR PGM'S





Fig. 1. Working curves for plane-grating monochromators with 300 lines/mm: (A) grating rotation at constant deviation  $|\alpha| + |\beta| = 2\theta = \text{constant}$  in inside diffraction order; (A') the same but outside diffraction order; (B) fix-focus mode,  $c_H = 2.25$ ; (B') fix-focus mode,  $c_H = 1/2.25$ .



# ACCESSIBLE RANGES OF THE SCAN PARAMETERS



- 300/mm grating
- $79^\circ < \theta < 88.4^\circ$
- $-89^\circ < \beta < -76^\circ$



# CONTRIBUTIONS TO THE SX700 RESOLUTION



Follath 1997, 2001

$$\frac{m\lambda}{d} = \sin\alpha + \sin\beta \quad \text{so} \quad \frac{\partial\lambda}{\partial\alpha} \frac{1}{\beta} = \frac{d \cos\alpha}{m}$$

$$\frac{d\lambda}{dq} = \frac{d\lambda}{d\beta} \frac{d\beta}{dq} = \frac{d \cos\alpha}{mr}$$

Fractional bandwidth due to the entrance slit is

$$\frac{\Delta\lambda}{\lambda}_{\text{source}} = \frac{\Delta q}{r} \frac{\cos\alpha}{m\lambda/d} = \frac{s}{r} A(c_{\text{ff}})$$

Angular error
Normalized angular dispersion

In the small - angle approximation

$$A(c_{\text{ff}}) = \sqrt{\frac{d}{m\lambda} \frac{2}{(c_{\text{ff}}^2 - 1)}}$$

NB - the grating changes the angular spread by  $1/c_{\text{ff}}$  ( $c_{\text{ff}}$  is the magnification)

$$\frac{\Delta\lambda}{\lambda}_{\text{source}} = \frac{s}{r} A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{exit slit}} = \frac{s}{r} c_{\text{ff}} A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{aberr}} = \frac{\Delta y}{r} c_{\text{ff}} A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{preopt}} = \sigma_{\text{po}} A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{focopt}} = \sigma_{\text{fo}} c_{\text{ff}} A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{grating}} = \sigma_{\text{gr}} (1 + c_{\text{ff}}) A(c_{\text{ff}})$$

$$\frac{\Delta\lambda}{\lambda}_{\text{total}} = \sqrt{\sum_i \left( \frac{\Delta\lambda}{\lambda}_i \right)^2} = \frac{1}{R}$$



# RESOLUTION CONTRIBUTIONS FOR A BESSY PGM



Follath 2001

BESSY U125/1-PGM

Foc optic 0.05 sec

Slit 10  $\mu\text{m}$

Grating 0.2 sec

Preoptic 0.25 sec

Source 50  $\mu\text{m}$

BESSY U125/1-PGM


Preoptic 0.25 sec

Source 50  $\mu\text{m}$

Foc optic 0.05 sec

Grating 0.2 sec

Slit 10  $\mu\text{m}$






# WORKING CURVES FOR $c_{ff}=2.25$ ETC



Petersen 1986





# BESSY CLSX700'S: SUMMARY



Table 1  
Operational collimated plane grating monochromators at BESSY II

Beamline	Energy range (eV)	Resolving power	Dedication
U125/1-PGM	15–600	100.000 <sup>a</sup>	low energy beamline
UE56/1-PGM	60–1500	100.000 <sup>a</sup>	circularly polarised, two chopped beams
UE56/2-PGM1	..	110.000 <sup>a</sup>	..
UE56/2-PGM2	..	110.000 <sup>a</sup>	..
U41-PGM	170–2000	6.000 <sup>b</sup>	high energy, high flux, small focus
U49/1-PGM	130–2000	10.000 <sup>b</sup>	PEEM
U49/2-PGM2	80–2000	10.000 <sup>b</sup>	PEEM
U180-PGM	30–1900 <sup>c</sup>	30.000 <sup>a</sup>	reflectometry, high spectral purity

<sup>a</sup> Aligned and determined with helium ( $2 - 1_3$ )-profile.

<sup>b</sup> Aligned and determined with  $N_2 : 1s \rightarrow \pi^*$  core hole absorption.

<sup>c</sup> Wiggler mode of the insertion device.




# BESSY U125/1-PGM BEAM LINE




Follath 2002

top view




Effective total length = 27 m

side view



distance between elements [mm]




distance to source point [mm]



# BESSY PGM RESOLUTION DEMONSTRATION




Picture from R. Follath, BESSY




- Doppler width 0.4 meV
- Monochromator contribution 0.65 meV
- Resolving power= $1.0 \times 10^5$
- Rotation increments
  - grating: 17 nrad
  - mirror: 9 nrad
- Measured with 1200/mm grating at  $c_{ff}=10$  to 12



# BESSY PGM'S: RESOLUTION SUMMARY







Parameters for this study:

150 lines/mm  
(low dispersion)

25nm groove depth

(operating from 75 eV to 800eV)

1200 lines/mm  
(high dispersion)

6nm groove depth

(operating from 300 eV to 2000eV)



# OPTICAL PATH FUNCTION: VARIED-LINE-SPACING GRATINGS



Now the groove spacing  $d(w)$  and number  $n(w)$  are functions of  $w$  and are also expressed as power series

$$d(w) = d_0 \left( 1 + v_1 w + v_2 w^2 + \dots \right) \quad n(w) = \sum_{i=1} n_{i00} w^i$$

So the optical path function becomes

$$F = \sum_{ijk} C_{ijk}(\alpha, r) w^i l^j z^k + \sum_{ijk} C_{ijk}(\beta, r) w^i l^j z^k + \frac{m\lambda}{d_0} n_{ijk}$$

$$n_{100} = 1$$

$$n_{200} = -v_1/2$$

$$n_{300} = (v_1^2 - v_2)/3$$

$$n_{400} = (-v_1^3 + 2v_1v_2 - v_3)/4$$

$$n_{ijk} = 0 \quad \text{if } i \text{ or } k = 0$$

Evidently the use of VLS can benefit aberrations of the form  $i00$ :

Defocus, coma, spherical aberration but not ones with  $j, k \neq 0$



# WHY DO VLS GRATINGS WORK?



The spherical - grating focus condition is now:

$$(F)_{200} = \frac{1}{2} w^2 \left( \frac{\cos^2 \alpha}{r} - \frac{\cos \alpha}{R} + \frac{\cos^2 \beta}{r} - \frac{\cos \beta}{R} - \frac{v_1 m \lambda}{d_0} \right) = 0$$

So setting  $R = \infty$ , using the grating equation and switching to grazing angles  $a, b$

$$\frac{1}{r} = \frac{\cos a - \cos b}{\sin^2 b} - \frac{1}{r} \frac{\sin^2 a}{\sin^2 b} = \frac{1}{r} v_1 \frac{-2 \sin \frac{a+b}{2} \sin \frac{a-b}{2}}{\sin^2 b} - \frac{\sin^2 a}{\sin^2 b}$$


Now approximating  $\sin a \approx a$  etc

$$\frac{1}{r} = \frac{1}{r} - v_1 \frac{a^2 - b^2}{2b^2} - \frac{a^2}{b^2}$$

Evidently if  $v_1 = \frac{-2}{r}$  then  $r = -r$



IDENTICALLY (for all values of  $a, b$  and thus  $\lambda$ . Moreover if, in addition,

$v_2 = \frac{1}{r^2}$  then coma and spherical aberration are corrected as well under the same small angle assumption






# DOES $\nu_1 = -r/2$ , $r' = -r$ GIVE OPTIMUM FOCUSING?



- No - the usual procedure is to write the focus equation twice for two wavelengths of exact focus and solve for  $r'$  and  $\nu_1$  - (this requires a fixed included angle or at least a fixed focusing principle) - the nominal focal distance then changes by a small amount (2% in the above case) - this gives locally better focus although not much improvement over a large range
- Note that all of the design studies in the literature except one (Koike 1995) consider only *fixed-included-angle* cases and offer, sometimes very good, although still approximate solutions but they have SGM-like disadvantages.
- Now that the variable-included angle schemes are available we have *exact* focus solutions such as CLSX700 or classical SX700 although the latter does not give control of the included angle to the user.
- Therefore retrofitting VLS to a standard SX700 makes some sense to give user control of the included angle however note that the focal length or the focal distance of the mirror would have to change.
- Conclusion: VLS role in general-purpose monochromators is becoming minor but it is still strong in spectrographs.



# VLS RETROFIT TO AN SX700




- VLS solution with  $c_{ff}$  chosen for maximum efficiency (between 1.5 and 2.5)
- Wavelengths of exact focus: 20 Å and 80 Å
- Source: size 50 μm and distance 15 m
- Beam height: 2mm
- Inside (+ve) order
- 1200/mm



# MES PROJECT CALCULATIONS




first order diffraction efficiency




First, second and third order efficiency calculations for a laminar grating 150/mm

second order diffraction efficiency



third order diffraction efficiency






# MES PROJECT CALCULATIONS



Source-size-limited resolving power for  
a source of size FWHM =  $40\text{ }\mu\text{m}$  at  
distance 12 m with a 150/mm grating






# MES PROJECT CALCULATIONS



cyl mirror rms slope error, R=7500 microrad

Cylindrical mirror slope error  
needed to get a resolving power of  
7500 with the 150/mm grating






# MES PROJECT CALCULATIONS





Grating slope error needed to  
get 7500 resolving power with  
150/mm grating







# Mirror cooling for high heat load at low energy (75eV)




The rms slope error ( $\mu\text{rad}$ ) of the pre-mirror corresponding to a resolving power  $R=7500$  (FWHM) from the 150l/mm grating.






# JENOPTIC SX700 MONOCHROMATOR






## New implementation of SX700 monochromator





# Monochromator chamber, vacuum test Oct 2001



M. R. Howells, Advanced Light Source



## monochromator scanning mechanism with integral cooling

